


Plan: Fusion^Prostate /PROSTATE
Structure Set: Fusion^Prostate /PROSTATE
Dose (Reference): Fusion^Prostate /PROSTATE/RTDose
Dose (Comparison): Fusion^Prostate /PROSTATE/RTDose
Report Generated: 2026/03/04 13:51:33
Reviewed by: _____

Dose Difference View

% Difference: 3% / Global Normalization (83.43261 Gy) / Lower TH = 10%
Distance Criterion: 2.0 mm Method: Gamma



Dose Grid Info

REFERENCE (PRIMARY) DOSE...

Dose Grid Resolution & Size (IEC Couch Coordinate System*):

Grid Spacing (X):	3.000 mm	X Range (mm):	[-187.500-->184.500]	*X along width of couch
Grid Spacing (Y):	3.000 mm	Y Range (mm):	[-076.500-->154.500]	*Y along long axis of couch
Grid Spacing (Z):	3.000 mm	Z Range (mm):	[-118.500-->133.500]	*Z normal to couch surface
Dose Min --> Max:	0.0000 --> 83.4326 Gy			


COMPARISON (SECONDARY) DOSE...

Dose Grid Resolution & Size (IEC Couch Coordinate System*):

Grid Spacing (X):	3.000 mm	X Range (mm):	[-187.500-->184.500]	*X along width of couch
Grid Spacing (Y):	3.000 mm	Y Range (mm):	[-076.500-->154.500]	*Y along long axis of couch
Grid Spacing (Z):	3.000 mm	Z Range (mm):	[-118.500-->133.500]	*Z normal to couch surface
Dose Min --> Max:	0.0000 --> 82.1108 Gy			

Plan: Fusion^Prostate /PROSTATE
Structure Set: Fusion^Prostate /PROSTATE
Dose (Reference): Fusion^Prostate /PROSTATE/RTDose
Dose (Comparison): Fusion^Prostate /PROSTATE/RTDose
Report Generated: 2026/03/04 13:51:33
Reviewed by: _____

Color	ROI Name	Volume (cc)	Quick Stat	Unit	Reference	Comparison	Diff (Abs / %)
	ref1 (Marker)	0.000	Point Dose**	Gy	59.999	60.004	0.005 Gy / 0.008%



Plan: Fusion^Prostate /PROSTATE
Structure Set: Fusion^Prostate /PROSTATE
Dose (Reference): Fusion^Prostate /PROSTATE/RTDose
Dose (Comparison): Fusion^Prostate /PROSTATE/RTDose
Report Generated: 2026/03/04 13:51:33
Reviewed by: _____

Comparison Parameters

% Difference: 3% / Global Normalization (83.43261 Gy) / Lower TH = 10%
DTA Criterion: 2.0 mm
Method: Gamma

Comparison Results: Main

Results: 84.7% Matching Rate (48 291 of 57 026 voxels, 6 793 high, 1 942 low)

Comparison Results: Per Anatomy ROI

ROI Name	Matching Rate (%)	Matching	High	Low	Total
patient	84.6	47875	6793	1942	56610
Bladder	92.6	1684	134	0	1818
GTV	79.6	1224	280	34	1538
dummy	80.3	10585	2357	233	13175
ref1		0	0	0	0

Plan: Fusion^Prostate /PROSTATE
Structure Set: Fusion^Prostate /PROSTATE
Dose (Reference): Fusion^Prostate /PROSTATE/RTDose
Dose (Comparison): Fusion^Prostate /PROSTATE/RTDose
Report Generated: 2026/03/04 13:51:33
Reviewed by: _____

----- Reference Data Sources -----

RTPlan

File: RTPLAN_2.16.840.1.114337.1.11224.1772428288.0.dcm

Dir: ~~C:\ProgramData\Sun Nuclear Patient Plans~~ PROSTATE_CCC

RTStructureSet

File: RTSTRUCT_2.16.840.1.114337.1.11224.1772428287.0.dcm

Dir: ~~C:\ProgramData\Sun Nuclear Patient Plans~~ PROSTATE_CCC

RTDose

File: RTDOSE_2.16.840.1.114337.1.11224.1772428288.1.CCC.dcm

Dir: ~~C:\ProgramData\Sun Nuclear Patient Plans~~ PROSTATE_CCC

Dose Grid Resolution & Size (IEC Couch Coordinate System*):

Grid Spacing (X):	3.000 mm	X Range (mm):	[-187.500-->184.500]	*X along width of
couch				
Grid Spacing (Y):	3.000 mm	Y Range (mm):	[-076.500-->154.500]	*Y along long axis of
couch				
Grid Spacing (Z):	3.000 mm	Z Range (mm):	[-118.500-->133.500]	*Z normal to couch
surface				
Dose Min --> Max:	0.0000 --> 83.4326 Gy			

----- Comparison Data Sources -----

RTDose

File: RTDOSE_2.16.840.1.114337.1.14324.1772511250.1.MC.dcm

Dir: ~~C:\ProgramData\Sun Nuclear Patient Plans~~ PROSTATE_CCC

Dose Grid Resolution & Size (IEC Couch Coordinate System*):

Grid Spacing (X):	3.000 mm	X Range (mm):	[-187.500-->184.500]	*X along width of
couch				
Grid Spacing (Y):	3.000 mm	Y Range (mm):	[-076.500-->154.500]	*Y along long axis of
couch				
Grid Spacing (Z):	3.000 mm	Z Range (mm):	[-118.500-->133.500]	*Z normal to couch
surface				
Dose Min --> Max:	0.0000 --> 82.1108 Gy			